

## LH HARMONY PRESO MASSAGE

### INTRODUCTION

Capable of simulating multiple massage techniques, the **LH HARMONY PRESO MASSAGE** bed is specially designed to try to relieve muscular discomfort.

Thanks to its wide range of possibilities it allows you to apply heat, massage or stimulate by acupressure the areas of the body most affected.

Another important advance of the device with respect to other models is its integrated computer processor. Combined with Chinese acupressure meridians, it will allow you to perform rhythmic stretching and body alloy movements on an air mattress.

It should be noted that **LH HARMONY PRESO MASSAGE** is not a medical device. Its use is recommended for all types of people, especially those who have physical and mental exhaustion, oppression, muscle tightness or any other type of discomfort.

Your body needs the best care and at Lufthous® we work so that you can enjoy the best possible resources from your own home.

**What are you waiting to find out?**

### OPERATION AND CHARACTERISTICS

The characteristics of the **LH HARMONY PRESO MASSAGE** bed make this product an innovative and practical resource.

Given the simplicity of its operation through its remote control, this device is extremely intuitive and easy to operate.

The controls give you access to four main functions;

- **Stretching:** With the right intensity to stimulate and prepare the body for subsequent exercises.

- **Massage:** After a long day of stress or physical exhaustion, there is nothing better than a massage to relax and relieve muscle tension.

- **Acupressure:** Simulating the famous Chinese technique, you will feel the fingers of an expert in each application.

- **Stretching of the spine:** An ideal function to help alleviate possible contraction and discomfort in the lower back.

Another point to highlight about **LH HARMONY PRESO MASSAGE** is its easy portability. Thanks to its light weight and side handles, you can take it anywhere!





## SPECIALLY RECOMMENDED FOR YOU

As mentioned above, **HARMONY PRESSO MASSAGE** is not a medical device. Its use is highly recommended, as long as the "Safety Warnings" are respected.

In the sporting field, this device can be a great ally. Perfect both before and after a workout.

Thanks to the **HARMONY PRESSO MASSAGE**, it is possible to alleviate the signs of tiredness. It also helps to improve internal sensations.

As for its application in relation to the aesthetic field, **LH HARMONY PRESSO MASSAGE** helps strengthen the elasticity and drainage of fat cells.

## LUFTHOUS® QUALITY

After years of work and effort, we have proven to be a benchmark of innovation, commitment and quality.

Each of our devices is the result of numerous investigations and conscientious field work. Only then, we can guarantee that you enjoy fabulous benefits and advantages in your own home that we could only dream of before.

All this wouldn't be possible without your trust. Thanks to you, our project grows and moves in the right direction.

What's our main objective? Your happiness. For this reason we thank you once again for the support and trust you have placed in each of our projects.

To a future together. To a better future.

## TECHNICAL SPECIFICATIONS

**Model:** HARMONY

**Product name:** LH PRESSO MASSAGE

**Voltage:** 100–240 V - 50/60 Hz  
12 VCC - 2 A

**Consumption::** 24 W

**Product size:** 120 cm × 60 cm × 9 cm (l × w × h)

**Product weight:** 4 kg



C/ Margarita 34, 36 y 38  
(POLIG. IND. EL LOMO)  
Humanes de Madrid. 28970 MADRID  
Teléfono: 902 040 252 / Fax: 91 279 15 19  
[902040252@lufthous.es](mailto:902040252@lufthous.es)  
<http://www.lufthous.es>

